



pH Tester for Skin

The health and beauty industry has been searching for a professional product that can determine the pH value of scalp, hair strands and skin quickly, accurately and inexpensively.

HANNA instruments® offers the *Skincheck*™, the first pH meter that makes it possible for you to easily measure pH in this application.

Skincheck™ measures pH in just a few seconds with a special flat-tip electrode specifically designed to measure pH on body surfaces. *Skincheck*™ has a measurement range that covers the full spectrum of hair and skin pH levels. Measurements are quickly displayed on a large easy-to-read LCD with an accuracy of 0.01.

Skincheck™ is an absolute must for Hair Care Specialists, Beauty Salons, Cosmeticians, Beauticians, Dermatologists & The Cosmetic Industry.



Specifications

| | HI 98109 (<i>Skincheck</i> ™) | HI 98110 (<i>Skincheck</i> ™1) |
|-------------------------|---|---------------------------------|
| Range | 0.00 to 14.00 pH | |
| Resolution | 0.01 pH | |
| Accuracy (@20°C/68°F) | ±0.2 pH | |
| Calibration | manual, 2 points | |
| Electrode | HI 1413S/50 (included) | HI 1413B/50 (included) |
| Battery Type / Life | 2 x 1.5V / approx. 3000 hours of continuous use | |
| Environment | 0 to 50°C (32 to 122°C); RH max 95% | |
| Dimensions | 66 x 50 x 25 mm (2.6 x 2.0 x 1.0") | |
| Weight (with electrode) | 70 g (2.5 oz.) | 90 g (3.2 oz.) |

Accessories

| | | | |
|-------------|--|-----------|---|
| HI 1413B/50 | Spare pH electrode for <i>Skincheck</i> ™1 with BNC connector and 1 m (3.3') cable | HI 7004M | pH 4.01 buffer solution, 230 mL bottle |
| HI 1413S/50 | Spare pH electrode for <i>Skincheck</i> ™ with screw connector | HI 7007M | pH 7.01 buffer solution, 230 mL bottle |
| HI 70004P | pH 4.01 buffer solution, 20 mL sachet, 25 pcs. | HI 7010M | pH 10.01 buffer solution, 230 mL bottle |
| HI 70007P | pH 7.01 buffer solution, 20 mL sachet, 25 pcs. | HI 70300M | Electrode storage solution, 230 mL bottle |
| HI 70010P | pH 10.01 buffer solution, 20 mL sachet, 25 pcs. | HI 7061M | Cleaning solution, 230 mL bottle |

Ordering Information

HI 98109 (*Skincheck*™) is supplied complete with HI 1413S/50 pH electrode, calibration screwdriver, batteries and instructions.

HI 98110 (*Skincheck*™1) is supplied complete with HI 1413B/50 pH electrode with 1 m (3.3') cable, calibration screwdriver, batteries and instructions.

For a complete range of calibration, cleaning and maintenance solutions, see section F.