


**HI 3618**  Connector

**HI 3618**  7-pin DIN



**Reference**

Single, Ag/AgCl

**Junction**

Ceramic, single

**Electrolyte**

KCl 3.5M + AgCl

**Max Pressure**

0.1 bar

**Range**

ORP T: -5 to 100°C (23 to 212°)

**Tip**

Platinum pin

**Temperature Sensor**

Yes

**Amplifier**

Yes

**Body Material**

Glass

**Cable**


5-pole; 1 m (3.3')

**Recommended Use**

Laboratory

**HI 3620**  Connector

**HI 3620**  7-pin DIN



**Reference**

Single, Ag/AgCl

**Junction**

Ceramic, single

**Electrolyte**

Gel

**Max Pressure**

2 bar

**Range**

ORP T: 0 to 80°C (32 to 176°F)

**Tip**

Platinum pin

**Temperature Sensor**

No

**Amplifier**

Yes

**Body Material**

Ultem®

**Cable**


7-pole; 1 m (3.3')

**Recommended Use**

Field applications, swimming pools

**HI 4619**  Connector

**HI 4619**  7-pin DIN



**Reference**

Single, Ag/AgCl

**Junction**

Ceramic, single

**Electrolyte**

Gel

**Max Pressure**

2 bar

**Range**

ORP T: -5 to 100°C (23 to 212°)

**Tip**

Gold pin

**Temperature Sensor**

Yes

**Amplifier**

Yes

**Body Material**

Glass

**Cable**

5-pole; 1 m (3.3')

**Recommended Use**

Strong oxidants, disinfection with ozone