

Reagents for Chemical Test Kits

| Code | Description |
|---|--|
| Soil analysis, Ned (N), ascorbic acid (P), tetraphenylborate (K), pH indicators (pH) | |
| Range: | Nitrogen (N): trace, low, medium, high Phosphorus (P): trace, low, medium, high Potassium (K): trace, low, medium, high pH: 4 to 9 pH |
| Reagents | |
| HI 3895-010 | Reagent kit for 10 tests (N, P, K, pH) |
| HI 3896-025 | Reagent kit for 25 tests (N, P, K, pH) |
| Sulphate, barium chloride method | |
| Range: | 20 to 100 mg/L |
| Reagents | |
| HI 38000-10 | Reagent kit for 100 tests (SO ₄ ²⁻) |
| Sulphate LR/HR, barium chloride method | |
| Range: | 100 to 1000 mg/L 1000 to 10000 mg/L |
| Reagents | |
| HI 38001-10 | Reagent kit for 200 tests (SO ₄ ²⁻) |
| Sulphite, iodometric titration method | |
| Range: | 0 to 20 mg/L 0 to 200 mg/L |
| Reagents | |
| HI 3822-100 | Reagent kit for 100 tests (Na ₂ SO ₃) |
| Zinc, extended range, zincon method | |
| Range: | 0.0 to 4.0 mg/L 0.0 to 20.0 mg/L |
| Reagents | |
| HI 38076-100 | Reagent kit for 100 tests (Zn) |
| Zinc, zincon method | |
| Range: | 0.0 to 3.0 mg/L |
| Reagents | |
| HI 3854-100 | Reagent kit for 100 tests (Zn) |

Accessories for Analysis

| Accessories | |
|-------------|--|
| HI 3824-99 | Ammonia color cube |
| HI 3830-99 | Bromine color cube |
| HI 3831-99 | Chlorine color cube |
| HI 3832-99 | Iodine color cube |
| HI 3833-99 | Phosphate color cube |
| HI 3834-99 | Iron color cube |
| HI 731340 | 200 µL automatic pipet |
| HI 731341 | 1000 µL automatic pipet |
| HI 731342 | 2000 µL automatic pipet |
| HI 731350 | Tip for 200 µL automatic pipet (25 pcs) |
| HI 731351 | Tip for 1000 µL automatic pipet (25 pcs) |
| HI 731352 | Tip for 2000 µL automatic pipet (25 pcs) |
| HI 740032P | Cap for 20 mL beaker (10 pcs) |
| HI 740034P | Cap for 50 mL beaker (10 pcs) |
| HI 740036P | Plastic Beaker, 50 mL (10 pcs) |
| HI 740037P | Plastic Beaker, 20 mL (10 pcs) |
| HI 740041 | 60 mL glass BOD bottle |
| HI 740048P | Plastic bottle, 120 mL (10 pcs) |
| HI 740049P | Cap for 120 mL plastic bottle (10 pcs) |
| HI 740132P | Dosing spoon (10 pcs) |
| HI 740142P | Graduated syringe (10 pcs) |
| HI 740144P | Tip for graduated syringe (10 pcs) |
| HI 740196P | Plastic vials (5 pcs) |
| HI 740220 | 25 mL glass test tube with cap, (2 pcs) |
| Demi-02 | Bottle for 2 L of demineralized water |
| Demi-05 | Bottle for 5 L of demineralized water |
| Demi-10 | Bottle for 10 L of demineralized water |

Demi-02, Demi-05 and Demi-10 Bottles for Demineralized Water

